

AMBU BAG ADULT CHILD INFANT

PRODUCT CODE:

Ambu Bag Adult: **STABAAB-03**
 Ambu Bag Child: **STABCAB-04**
 Ambu Bag Infant: **STABIAB-05**

PRODUCT FEATURES:

- * Provide assisted ventilation to people who are either not breathing or are having trouble breathing.
- * The bag needs to be compressed to force a volume of air into the lungs.
- * Consists of a bag, adapter and one-way valve. The bag is the main part of the device.
- * It can be attached directly to an endotracheal tube or a mask using the adapter.
- * The bag also has a one-way valve which allows the oxygen to flow to the patient and prevents the patient from breathing in his exhaled air.
- * Some bags may have a reservoir, which increases the amount of oxygen the patient receives.



PRODUCT SPECIFICATION:

Specifications: Infant/Neonate Pediatric Adult
 Body weight: ≤ 10kg 10 kg~33 kg ≥33 kg
 Stroke volume: 150ml/100ml 350ml 800ml
 Resuscitator volume: 340ml 680ml 1630ml
 Resuscitator weight: 350g 410 g 600 g
 Pressure limiting valve: 40cmH2O 40cmH2O 60cmH2O

Dead space: 7ml 7ml 7ml
 Bag reservoir volume: 1600ml/600ml 1600ml 2000ml
 Inspiratory resistance: <5 cmH2O (at 50 L/min)
 Expiratory resistance: <5 cmH2O (at 50 L/min)
 Patient connector: ISO Ø22/15 mm
 Recommended Operating Temp.: -18°C ~ +50°C
 Storage: -40°C ~ + 60°C


Certification:



FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 www.abronn.com

 info@abronn.com

 +971(0) 48844860

 +971 508804579