

## AIR SPLINT

**PRODUCT CODE: STAISAB-02**

**PRODUCT SPECIFICATION:**

- \* Immobilize arm and leg fractures, with uniform pressure on the limb and to guarantee stability and comfort for the patient.
- \* It consists of an inflatable cylinder that can be closed at either end and becomes rigid when filled with air under pressure.
- \* It is easy to apply and remove.
- \* It has a push-pull valve for hassle-free inflation.
- \* Easy to use and economical.




**Certification:**



FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 [www.abronn.com](http://www.abronn.com)

 +971(0) 48844860

 [info@abronn.com](mailto:info@abronn.com)

 +971 508804579