

BLOOD PRESSURE MONITOR

PRODUCT CODE: STBPMAB-10

PRODUCT SPECIFICATION:

- * Display: Digital liquid crystal display
- * Accuracy: static: $\pm 3\text{mmHg} (\pm 0.4\text{kPa})$
- * Measurement mode: Oscillometric method
- * Measurement range: Pressure: $0\text{kPa}-37.3\text{kPa}$ ($0\text{mmHg}-280\text{mmHg}$)
- * Resolution: 0.1kPa (1mmHg)
- * Power supply: DC 6V (1.5V four AAA batteries)
- * Storage temperature and moisture: Temperature ($-10^{\circ}\text{C}-55^{\circ}\text{C}$)
Moisture: $\leq 10\%-85\%RH$)
- * Using temperature and moisture: Temperature ($5^{\circ}\text{C}-40^{\circ}\text{C}$)
Moisture: $\leq 15\%-85\%RH$)
- * Weight: each one is about 220 grams
- * Dimention: $126 \times 100 \times 53\text{mm}$ (not include Wristbands)
- * Power supply: AC/DC 6V US plug not include adapter/
Lithium battery




Certification:




FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 www.abronn.com

 +971(0) 48844860

 info@abronn.com

 +971 508804579