

CPR MASK

PRODUCT SPECIFICATION:

- * The CPR pocket mask is a non-invasive, latex free mask with inflatable cuff and an integral non-rebreathing valve and efficient filter to prevent direct mouth contact with the patient's face when providing mouth-to-mouth resuscitation.
- * It can be used on adult, child and infant.
- * Supplied in a compact plastic box with a corrugated tube for an easier use.




Certification:



FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 www.abronn.com

 +971(0) 48844860

 info@abronn.com

 +971 508804579