

EMOX POWDERED EMERGENCY OXYGEN

PRODUCT CODE: STEPEAB-76

PRODUCT SPECIFICATION:

- * It is a portable, reusable device that produces 99.7% pure O₂ through a chemical reaction of 450ml of tap water, 300g of sodium percarbonate and 7g of manganese catalyst in a polypropylene plastic cylinder.
- * The oxygen is created within 15 to 30 seconds, is humidified through a water chamber and lasts for approximately 15 minutes.
- * The powders have an indefinite shelf life if stored at room temperature.
- * Emox powdered oxygen is used to administer pure O₂ with a respirator face mask in emergency situations such as traumatic accidents, sporting & outdoor first aid and in medical application before receiving professional healthcare, or as part of oxygen therapy.




Certification:




FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 www.abronn.com

 +971(0) 48844860

 info@abronn.com

 +971 508804579