

EVACUATION CHAIR

PRODUCT CODE: STEVCAB-83

PRODUCT FEATURES:

- * With track, only one person can push down the stretcher.
- * It is made of high-strength powder-coated aluminum alloy, PVC seat.
- * Telescoping backrest
- * Foldable, easy for storage

PRODUCT SPECIFICATION:

High Position (LxWxH): 92x53x138cm
Low Position (LxWxH): 35x53x115cm
N.W.: 11kg G.W.: 13 kg
Loading capacity: 160kg




Certification:




FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 www.abronn.com

 info@abronn.com

 +971(0) 48844860

 +971 508804579