

## HEMOGLOBIN TEST KIT

**PRODUCT CODE: STHESAB-39**

**PRODUCT CONTENT:**

- 1 HB Meter
- 1 Owner's Manual
- 1 Storage Bag
- 3 AAA Batteries
- 1 Code Chip

**PRODUCT SPECIFICATION:**

Blood Sample: Capillary and venous whole blood  
 Sample Volume: 10uL  
 Measuring Time: Less than 15 seconds  
 Measuring Range: 4.5-25.6g/Dl  
 Scale: switchable in g/dL,g/L,mmol/L  
 LCD display: 37x39mm  
 Battery type: 3 AAA alkaline batteries  
 Unit dimension: 12.6x5.7x2.5cm  
 Unit weight: Approx.185grams  
 Packing: 1PCs/Giftbox; 24PCs/Cartron

\*Frequent testing helps to improve and control blood glucose levels.





**Certification:**



FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 [www.abronn.com](http://www.abronn.com)  
 [info@abronn.com](mailto:info@abronn.com)

 +971(0) 48844860  
 +971 508804579