

PORTABLE SUCHION MACHINE

PRODUCT CODE: STPSMAB-57

PRODUCT FEATURES:

- * It is designed with oil free pump for environmental protection and low noise.
- * Overflow protection and inside bottle make compact design, lightweight, small size and is convenient for transportation and safe use.

PRODUCT SPECIFICATION:

Voltage: AC220V±10%

Input power: 90VA

Noise: ≤65dB

Max negative pressure value: ≥0.075MPa

Negative pressure range: 0.02 MPa to Max negative pressure

Pumping rate: 15L/min

Jug Volume: 1000ML

Continuous Working Time: 0-30 min

Rest Time: ≥30min

Bottle volume: 1000 ml

Package: 2pcs/carton

Box Size: 54.5*37*32cm (2pcs)

Accessory: Soft tube(1pc, 2m long, 2mm thick, 12mm in diameter),
suction tube(2pcs, 1 for adults, 1for children, 55cm per pc),

fuse (2pcs, 1.5A & 3A), filters(2pcs), (Color box for retail purpose)

Gross Weight: 16 kg (2pcs/carton)





Certification:



FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 www.abronn.com

 +971(0) 48844860

 info@abronn.com

 +971 508804579