

PORTABLE SUCTION MACHINE

PRODUCT CODE: STPSMAB-57

PRODUCT SPECIFICATION:

Voltage: **AC220V±10%, 50Hz±2%, 12VDC**

Input power: **90VA**

Noise: **≤65dB**

Max negative pressure value: **≥0.075MPa**

Negative pressure range: **0.02 MPa to Max negative pressure**

Pumping rate: **15L/min**

Jug Volume: **1000ML**

Continuous Working Time: **0-30 min**

Rest Time: **≥30min**

Bottle volume: **1000 ml**

Package: **2pcs/carton**

Box Size: **54.5*37*32cm(2pcs)**

Accessory: Soft tube(1pc, 2m long, 2mm thick, 12mm in diameter), suction tube(2pcs, 1 for adults, 1for children, 55cm per pc), fuse (2pcs, 1.5A & 3A), filters(2pcs), charger(1pc, used in car)

Gross Weight: **16 kg (2pcs/carton)**




Certification:




FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 www.abronn.com

 +971(0) 48844860

 info@abronn.com

 +971 508804579