

## SCOOP STRETCHER


### **PRODUCT CODE: STSCSAB-62**

### **PRODUCT FEATURES:**

Product Size (L x W x H):  
Max. Size: 210 x 44 x 6cm  
Min. size: 120 x 44 x 6cm  
Packing size: 170 x 46x 8cm  
N.W.: 8.5kg G.W.: 10.2kg  
Load bearing: 159kg

### **PRODUCT SPECIFICATION:**

- \* Made of high-strength aluminum alloy materials.
- \* A separately-type stretcher to transfer fractured and serious patients.
- \* A hinged and joined arrangement, which is easy to set the stretcher into two parts, is set in the middle at both ends of the stretcher.
- \* Strengthen the patients at original place to reduce the further hurt of the patients.
- \* Put the patients in the stretcher quickly & easily.
- \* Can be taken out from the back of the patient without moving the patient.
- \* Adjust the length according to patients' body.
- \* A narrow-frame foot structure.
- \* Light-weighted, small-sized, easy-carried, use-safely and easy for sterilization and clean.




**Certification:**




FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 [www.abronn.com](http://www.abronn.com)

 [info@abronn.com](mailto:info@abronn.com)

 +971(0) 48844860

 +971 508804579