

## STAIR CHAIR

### **PRODUCT SPECIFICATION:**

- \* This stair stretcher is made of high-strength aluminum alloy materials. It is an ideal facility for emergency staff to transfer patient to go up and down stair.
- \* Four wheels enable the patients to be moved over most surfaces with minimal effort.
- \* The stretcher approves the patient safety during transport process because of its 2 pieces of restraint straps with quick-release buckle.
- \* Under the front of the stretcher, it has two flexible uplift rods.
- \* The structure can be folded and it saves the space in ambulance car.
- \* It characterized by its being light-weighted, small sized, easy-carried, use-safely and easy for sterilization and clean.

### **PRODUCT FEATURES:**

Technical Features:

Product Size: (L\*W\*H): 90x 51 x91cm

Folding Size: 90x17x59xcm

Packing Size: 96x16x54cm

NW: 8kg GW: 9.2kg

Load bearing: Less than 160kg




**Certification:**



FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 [www.abronn.com](http://www.abronn.com)

 [info@abronn.com](mailto:info@abronn.com)

 +971(0) 48844860

 +971 508804579