

SAM SPLINT

PRODUCT CODE:

Sam Splint (large): **STSSLAB-59**
 Large Size: **11cmx 100cm**
 Sam Splint (medium): **STSSMAB-60**
 Medium Size: **11cmx75cm**
 Sam splint (small): **STSSSAB-61**
 Small Size: **11cmx50cm**

- * Built from a thin core of soft aluminum alloy sandwiched between two layers of closed cell foam.
- * SAM Splint bends into any of three simple curves to give the splint strength to immobilize a fractured or injured limb.



Certification:



FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 www.abronn.com
 info@abronn.com

 +971(0) 48844860
 +971 508804579