

## SHORT SPINAL BOARD

**PRODUCT CODE: STSBSAB-65**

**PRODUCT FEATURES:**

- \* Spine board is commonly used for safely transporting seriously injured and immobile.
- \* It can float on water, allowed X-rays and also used together with head immobilizer.
- \* Made of high strength engineering plastic, this spine board is enduring, anti-aging, weight loading and easy for storage.
- \* It is a must have in first-aid equipment for medical centers.

**PRODUCT SPECIFICATION:**

Product size (L x W x H): **138 x 35 x 5cm**  
Packing size (1pc/carton): **140 x 38 x 7cm**  
N.W.: **4.6kg** G.W.: **6kg**  
Load bearing: **159kg**




**Certification:**




**FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE**

 [www.abronn.com](http://www.abronn.com)

 **+971(0) 48844860**

 [info@abronn.com](mailto:info@abronn.com)

 **+971 508804579**