

TRACTION SPLINT

PRODUCT FEATURES:

- * Relieve the tissue pressure on fractured leg avoiding further hemorrhages, shock and secondary nervous, vascular, bone and muscular damages.
- * Made of aluminum alloy.
- * Adjustable length, and nylon leg straps automatic traction mechanism.
- * Set of adult and pediatric sizes together.
- * With carry bag.

PRODUCT SPECIFICATION:

Product Features (ADULT):

- * Product size (LxWxH): 138x18x17cm
- * Folded size (LxWxH): 93x7x18cm
- * Packing size (1pc/carton): 95x8x22cm
- * N.W.: 2kg G.W.: 3kg

Product Features (CHILD):

- * Product size (LxWxH): 123x18x17cm
- * Folded size (LxWxH): 86x7x18cm
- * Packing size (1pc/carton): 88x8x22cm
- * N.W.: 1.8kg G.W.: 2.8kg




Certification:




FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 www.abronn.com

 info@abronn.com

 +971(0) 48844860

 +971 508804579