

VACUUM MATRESS w/ PUMP

PRODUCT FEATURES:

The removable strapping and handle system reduces the areas that could trap fluids. One valve evacuates the entire unit. Lightweight, durable immobilization system designed to fold down to a compact size for storage and sets up quickly for immediate use.

Consist of :

- * Mattress
- * Plastic vacuum pump
- * Nylon carry case

PRODUCT SPECIFICATION:

- * Immobilizes the patient in a secure position
- * Conforms to any body shape
- * The wrap-around straps with buckles to provide maximum patient security
- * Length: 202 cm, Width: 92 cm Height: 5 cm Weight: 7 kg
- * Load Capacity: 125kg




Certification:



FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 www.abronn.com

 info@abronn.com

 +971(0) 48844860

 +971 508804579