

VACUUM SPLINT

PRODUCT CODE: STVASAB-73

PRODUCT FEATURES:

Vacuum splints provide the best possible support and versatile extremity splinting system. The design is ideal for dislocated shoulders, suspected fracture of arms, legs, and joints, or infant immobilization. They suite a wide variety of emergency applications.

- * Insert the 5-in-1 connector securely into valve, requires only one hand. When connector is withdrawn, valve closes automatically.
- * To release vacuum after use, simply squeeze the valve.

PRODUCT SPECIFICATION:

- * Small splint
- * Medium splint
- * Large splint
- * Aluminum vacuum pump
- * Nylon carry case




Certification:




FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 www.abronn.com

 +971(0) 48844860

 info@abronn.com

 +971 508804579