

WALL MOUNTED SPHYGMOMANOMETER

PRODUCT CODE: STWMSAB-74

PRODUCT SPECIFICATION:

- * It has a large round face.
- * This model includes an adult cuff + 2-tube bladder and a long coil tube.
- * Measurement ranges from 0 mmHg up to 300 mmHg.
- * The inflation of cuff and air release are done manually.
- * Also can be used with a standard stethoscope.




Certification:




FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 www.abronn.com

 +971(0) 48844860

 info@abronn.com

 +971 508804579