

WHEEL CHAIR

PRODUCT FEATURES:

FRAME:	Double cross braces, liquid painting
FRONT WHEEL:	8"x1.25", PU tires
REAR WHEEL:	24", spoke wheel, PU tires
ARMREST:	Flip back slope armrest
FOOTREST:	Swing in & away footrest
SEAT AND BACKREST:	Nylon with sponge
FRONT FORK:	Plastic
FOOTPLATE:	Plastic
PACKAGE:	24"Rear Wheel: 119PCS/20"GP, 302pcs/40"HQ
	12"Rear Wheel: 152PCS/20"GP, 360pcs/40"HQ
OPTION:	Rear wheel: 12", PU tires; Height adjustable armrest
	Backrest half folding; attendant brakes
SPECIAL:	Quick release



SIZE	16"	18"
Overall Length (cm)	103	103
Overall Width (cm)	60	65
Overall Height (cm)	92	92
Overall Width Folded	30	30
Seat Width (cm)	43	46
Seat Depth (cm)	41	41
Back Height (cm)	43	43


Floor to seat Height	49	49
Weight Capacity (kg)	135	135
Front Caster Spec.	8"	8"
Rear Wheel Spec.	24"	24"
G./N. Weight (kg)	15.8/13.8	16.3/14.3
Carton Size (cm)	82x31x93	82x31x93

Certification:




FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 www.abronn.com

 +971(0) 48844860

 info@abronn.com

 +971 508804579