

AUTOMATIC LOADING STRETCHER

Light Weight

Durable

Safe Usage

Easy for
Sterilization

PRODUCT SPECIFICATION:

- * This is made of high-strength aluminum alloy thick-wall pipes.
- * It is supplied with foamed cushion and the adjustable back making the wounded feel comfortable.
- * This machine is very convenient; one person is enough to send the patient (with this stretcher) to the ambulance car.
- * The mechanical foldaway legs are controlled by the handles on the both ends of the stretcher, & the legs can be folded under the stretcher automatically when loading into the vehicle.
- * Also, the stretcher can be locked by the fixing devices when on the ambulance car.



PRODUCT FEATURES:

High Position: 185 x 55 x 92 cm

Low Position: 185 x 55 x 25 cm

Max Angle of the back: 75°

Distance between the front and rear wheel: 108cm

Maximum Load Bearing Capacity: 159kg

N.W.: 46kg G.W.:52kg

Certification:



FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

www.abronn.com

+971(0) 48844860

info@abronn.com

+971 508804579