

## FOLDABLE STRETCHER



**PRODUCT CODE: STFOSAB-30**

**PRODUCT FEATURES:**

- \* The stretcher is made of high strength aluminum alloy material and the waterproof leather fabric.
- \* Stretcher can be folded in length and width, easy for transportation and storage.
- \* The fabric is fire proof, water proof, anti-cracking and easy to clean.
- \* With safety belts, stirup feet and hinged spread bar.

**PRODUCT SPECIFICATIONS:**



Product Size (LxWxH): **208x 55x13cm**  
 Folding Size(LxWxH): **104x17x9cm**  
 Packing Size (2pcs): **108x18.5x18cm**  
 N.W.: **4.7kg** G.W. : **5kg**  
 Load Bearing: **159kg**

**Certification:**



FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 [www.abronn.com](http://www.abronn.com)  
 [info@abronn.com](mailto:info@abronn.com)

 +971(0) 48844860  
 +971 508804579