

HEAD IMMOBILIZER

PRODUCT CODE: STHEIAB-38

PRODUCT FEATURES:

- * Made of high-density plastic materials, large ear holes for monitoring the patient's ear canal
- * It can be used with scoop stretcher and spine board at the same time
- * Minimal interference with X-ray, MRI, or CT scanning procedures
- * Water-resistant plastic coating makes it easy to clean and prevents bacterial growth on it.

PRODUCT SPECIFICATION:

Padding : 40 cm x 25 cm x16cm
Two head blocks: on the side: 13.8cm
on the bottom: 9cm
on the top: 5.2 cm
External diameter of the hole: 6 cm
Packing size (L x W x H): 41*27*20cm
N.W: 1kg G.W: 1.5kg




Certification:




FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 www.abronn.com

 +971(0) 48844860

 info@abronn.com

 +971 508804579