

## LONG SPINAL BOARD

**PRODUCT CODE: STSBLAB-64**

**PRODUCT FEATURES:**

- \* This spine board with safe belts is made of PE material with no discharge contaminator, and firmness to wear.
- \* This plastic spine board provided by our company also can be used as a floatable device.
- \* The structure is allowed X-rays because of its completely
- \* translucent.  
It is compatible with most head immobilization devices and
- \* strap mechanisms.  
It is mainly used for hospitals, sports, ambulance and outdoors to carry patients and wounded person.

**PRODUCT SPECIFICATION:**

Product size (L x W x H): **183 x 45 x 5cm**  
Packing size (1pc/carton): **185 x 47 x 7cm**  
N.W.: **7.5kg** G.W.: **9kg**  
Load bearing: **159kg**




**Certification:**




FZS1 BG03-BG04, Jebel Ali Free Zone, P. O. Box 487449 Dubai, UAE

 [www.abronn.com](http://www.abronn.com)

 +971(0) 48844860

 [info@abronn.com](mailto:info@abronn.com)

 +971 508804579