

# MERCEDES SPRINTER



★ HIGH QUALITY    ★ ON-TIME DELIVERY    ★ COMPETITIVE PRICE



Plot No. S21504, Jebel Ali Free Zone South,  
P.O. Box 487449, Dubai, UAE



+971 4 884 4860  
+971 50 880 4579



sales@abronn.com  
www.abronn.com

# About Abronn FZE



Abronn FZE is one of the leading provider in healthcare solutions, medical units, special vehicles and accessories in the United Arab Emirates and around the world. Our company primarily offers wide range of ambulances, emergency transport, mobile clinics, field hospitals, special vehicle conversion and other mobile solutions. Additionally, we supply medical devices, hospital and emergency equipment. We also cater equipment for fitness, sports, beauty, leisure and outdoor.

We are recognised as CE certified manufacturer for our products from vehicles, machinery and medical devices. Abronn consistently delivers high-quality service tailored in the exact needs of our customers. Adding value to our products and services, Abronn FZE is accredited certified with ISO 9001: Quality Management System (QMS), ISO13485: Medical Devices, ISO14001: Environmental Management System, and ISO45001: Occupational Health & Safety.

► **Design & Technology Expertise**  
By combining our designers knowledge, flair and experience in 3D design, our design solution are innovative, practical and cost effective.

► **Product Availability**  
With dedicated warehouse for ambulance vehicles, medical equipment and other devices, our customers are able to benefit from availability, delivered fast when its needed, on time and in full.

► **Competitive Pricing**  
Our products are best in the market, we see to it that the cost is considerable providing best quality, safe and efficient solution for our customer.

► **After Sales Support**  
Our committed team works with our customer from the beginning, ensuring the products are in best condition when transported, support every concern until post delivery; sale/technical queries.

# EXTERIOR FEATURE




**LED LIGHT BAR**  
(WITH SIREN AND  
SPEAKER & MICROPHONE)



**FRONT STROBE LIGHTS**  
(COLOR OPTIONS)



**SIDE STROBE LIGHTS**  
(COLOR OPTIONS)



**AC UNIT W/ SEPARATE  
COMPRESSOR**



**STICKER SIGN**  
(REFLECTIVE)



## FRONT SIDE VIEW



① **LED LIGHT BAR**  
(WITH SIREN, SPEAKER  
& MICROPHONE)



② **REAR STROBE  
LIGHTS**



③ **STICKER SIGN**  
(REFLECTIVE)



## BACK SIDE VIEW



**INTERIOR FEATURE**

**A PATIENT COMPARTMENT**

- Insulation and covering of the patient compartment
- Partition with sliding window
- Cabinet for medical equipment
- Oxygen cylinders cabinet
- Flooring: marine plywood covered with vinyl
- Squad bench with under seat storage, backrest & seat belt
- Doctor seat

**B ELECTRICAL SYSTEM**

- Electric wiring with insulation
- Electric control panel/fuse box
- Roof type ventilation fan
- 2x Socket – 12V & 2x Outlet 220V
- Shore line (10m & 220V)
- Inverter
- LED lighting for the patient compartment
- Reverse alarm
- LED light bar with siren, speaker and microphone
- Strobe lights at the back of the vehicle
- LED Loading light
- WEBASTO, Air conditioning unit for the patient compartment, with separate compressor

**C OXYGEN SYSTEM**

- 2x 10L Oxygen cylinder, 1x oxygen regulator, outlets
- (3x 12V) (2x 220V) & 2x flow meter humidifier
- Suction pump, regulator, outlet & suction jar
- Oxygen system related hoses

**D PATIENT TRANSPORT EQUIPMENT**

- Automatic Loading Stretcher
- Stretcher Platform
- Foldable Stretcher
- Stair Chair Stretcher

**E MEDICAL EQUIPMENT**

- |                         |                         |
|-------------------------|-------------------------|
| ■ First Aid Kit         | Others Accessories:     |
| ■ Wall mounted          | ■ Fire Extinguisher     |
| ■ Sphygmomanometer      | ■ Handles               |
| ■ Stethoscope           | ■ Trash Bin             |
| ■ Syringes              | ■ Hammer                |
| ■ Manual Suction Device | ■ Side reflective lines |
| ■ Bandages              | ■ IV Hooks              |
| ■ Sam Splint            | ■ Sharps container      |
| ■ Emergency Blanket     | ■ Wash Basin (Optional) |